



WiZARD 310/320

Nasal Mask & Full Face Mask

- Contour fitting delivers better facial pressure redistribution
- The flexible frame without a forehead pad provides further comfort
- The snap on cushion allows instant assembly of mask and easy cleaning
- Special air guiding vents create a silent sleep environment by minimizing undesired airflow noise

WiZARD 310 & 320 Nasal Mask & Full Face Mask

The WiZARD 310 & 320 have a flexible frame material and contour design that bring further comfort and to help you sleep through the night



WiZARD 310
Nasal Mask



WiZARD 320
Full Face Mask



► Soft Forehead Support

Reduces the chance of losing the forehead pad and the nuisance of assembly and cleaning

► Contour Fitting Design

The ergonomic cushion with multiple-pressure-redistribution design provides better facial support to each individual



► Quick Snap on Design

The easy snap on cushion design makes quick mask replacement possible



► Advanced Air Guiding Vent

The oriented exhalation flow design prevents air from blowing onto the sleeping partner and minimizes high frequency noise



Soft Forehead Support
The soft and flexible forehead support structure with less skin contact provides further comfort and retains stability

Flexible Mask Frame
The flexible and lightweight frame offers better fitting and great seal

Breathable and Durable Headgear
Ensures a comfortable and effective OSA treatment

Safety Valve Structure
Designed to avoid asphyxia in case of CPAP machine power failure (for WIZARD 320 only)

Soft Connection Tubing
Offers better flexibility for freedom of movement during sleep

Quick Release Tubing
An easy one step removal design for convenient tubing connection

Headgear Clip
One buckle operation for simple mask fitting

Ordering Information (For North America Only)

Order No.	Description	Order No.	Accessories
SM07001	WIZARD 310 Nasal Mask with Headgear, Large	SM00042	WIZARD 310/320 Buckles
SM07002	WIZARD 310 Nasal Mask with Headgear, Medium	SM00043	WIZARD 310/320 Headgear with Buckles
SM07003	WIZARD 310 Nasal Mask with Headgear, Small	SM00044	WIZARD 310 Nasal Mask Cushion, Large
SM08001	WIZARD 320 Full Face Mask with Headgear, Large	SM00045	WIZARD 310 Nasal Mask Cushion, Medium
SM08002	WIZARD 320 Full Face Mask with Headgear, Medium	SM00046	WIZARD 310 Nasal Mask Cushion, Small
SM08003	WIZARD 320 Full Face Mask with Headgear, Small	SM00047	WIZARD 320 Full Face Mask Cushion, Large
Order No.	Accessories	SM00048	WIZARD 320 Full Face Mask Cushion, Medium
SM00040	WIZARD 320 Elbow Connector with frame	SM00049	WIZARD 320 Full Face Mask Cushion, Small
SM00041	WIZARD 310 Elbow Connector with frame	SM00016	Small Silicon Tubing

HEADQUARTERS

Wellell Inc.

No. 9, Minsheng St., Tucheng Dist.,
New Taipei City 236044, Taiwan
T +886 2 2268 5568
F +886 2 2268 6525
marketing@wellell.com

PAN-ASIA

TAIWAN

Wellell Inc.
No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 236044, Taiwan
T +886 2 2268 5600 / 0800 054 133
F +886 2 2268 6526
customer@wellell.com

STURDY INDUSTRIAL

Sturdy Industrial Co., Ltd.
168, Sec.1, Zhongxing Road, Wugu
District, New Taipei City 24872, Taiwan
T +886 2 2981 5245
F +886 2 2984 8208
foreign@sturdy.com.tw
www.sturdy.com.tw

CHINA

Wellell (Kunshan)Co., Ltd.
No. 1368, ZiZhu Road, Kunshan,
Jiang Su Province, China
T +86 512 8617 1660 Ext 1150
F +86 512 8617 8663
customer.cn@wellell.com

THAILAND

Apex Medical (Thailand)
No.111/152 Soi Pho Kaoe Yak 19,
Nawamin, Beung Kum,
Bangkok 10240,
Thailand.
T +66 02 018 2140
sales.thai@wellell.com

PAN-EUROPE

SPAIN

Apex Medical S.L.
Elcano 9, 6ª planta 48008
Bilbao. Vizcaya. Spain
T +34 94 4706408
F +34 94 4706409
info@apexmedical.es

UK

Wellell UK Limited
Unit 33, Great Western Business Park
McKenzie Way, Worcester ,WR4 9GN
United Kingdom
T +44 (0) 190 5774 695
sales@wellell.com

FRANCE

Wellell France S.A.S
4 Boulevard de la Chanterie
49124 Saint-Barthelemy d' Anjou
France
T +33 2 72 79 27 13
F +33 2 85 73 07 67
info.france@wellell.com

GERMANY

SLK Vertriebsgesellschaft mbH
Am Herdicksbach 18, 45731 Waltrop
Germany
T +49 (0) 231-925360-0
F +49 (0) 231-925360-29
info@slk-gmbh.de
www.slk-gmbh.de

AMERICAS

USA

Wellell America Corp.
927 Mariner Street,
Brea, CA 92821, USA
T +1 714 671 3818
F +1 714 494 8173
sales.usa@wellell.com

Distributed by